

4 TRIZONE

**Disposable Endoscopic Articulating Reload With
FUNNEL Technology**

SURGICAL STAPLERS

Laposcopic stapler

TRIZONETM

Disposable Endoscopic Articulating Reload With FUNNEL Technology



Large opening at the front of the reloads, convenient for organization.



The "drop-shaped" nail groove ensures that the B-shaped staples are formed into nails which improves the clinical suture molding effect.



360° rotation and 45° left and right swings can reach the target site more conveniently and accurately.



The reinforced fixed nail anvil improves the tissue closing force and ensures an effective cutting and closing length. It is easier to pass through a narrow area and is more suitable for delicate surgery.



The double arrow indicates the installation position. The assembly and disassembly of reloads are simple and convenient, and the surgical blade is pulled back smoothly.



The brand-new guiding technology of the reload makes the tissue suture smooth. At the same time, the safety design of the empty reload prevents the used order from being accidentally fired.



TRIZONE™


Disposable Endoscopic Articulating Reload With FUNNEL Technology

FEATURES

- Gradual compression for vascular, medium & thick tissues.
- Better blood supply, Better tissue thickness adaptability, Less tissue crush injury.
- The design of six rows cross staplers keeps the tissue safe.
- Single use blade, Changed with reload, Ensure the blade is sharp and reduces tissue damage.
- Equipped with the safety inter lock device to prevent the used reload from being triggered again. Thus ensures safety in surgical operation.
- Isometric height retainer structure, Made the staple mold consistently.
- Different altitude for staple, The anvil is suitable for pressing various tissue thickness, effectively clamping and firing.

SPECIFICATION

PRODUCT CODE	HEIGHT (MM)	LENGTH(MM)	NO OF STAPLES
FZHPL-45W	2.4/2.6/3.0	46	66
FZHPL-45P	3.4/3.6/4.0		
FZHPL-45B	4.4/4.6/4.8		
FZHPL-60W	2.4/2.6/3.0	61	90
FZHPL-60P	3.4/3.6/4.0		
FZHPL-60B	4.4/4.6/4.8		



SCAN THIS CODE
TO VISIT OUR WEBSITE

Address:-A-205 & 208,Premium House,Nr Gandhigram,
Railway Station,Ashram Road,Ahmedabad-380006.

☎ info@fourzone.in | sales@fourzone.in