

C-ANAZONETM


Disposable Circular Stapler with MAC Technology

SURGICAL STAPLERS

Open Stapler

C-ANAZONETM

Disposable Circular Stapler with MAC Technology



C-ANAZONE™

Disposable Circular Stapler with MAC Technology

FEATURES


- MAC Technology - Minimum gap between anvil and cartridge to avoid the damage to tissue caused by over pressure.
- Integrally shaping anvil to ensure the shape of staples.
- Large diameter ratio, Avoids anastomotic stenosis effectively.
- Wider viewport for precise examination.
- Reduce post-operative hemorrhage with appropriate arrangement of staples.
- Super-thin anvil, More convenient for procedure.
- Adjustable thickness of tissue, Meet all kinds of need.

SPECIFICATIONS

PRODUCT CODE	DIAMETER(MM)	CUTTING INNER DIAMETER (MM)	NO OF STAPLES
FZCIR-21	21.8	13.2	18
FZCIR-23.5	24.3	15.2	18
FZCIR-25.5	26.2	17.2	20
FZCIR-28.5	29.2	20.2	24
FZCIR-31.5	32	22.2	30
FZCIR-33.5	33.8	24.2	32

APPLICATION

- It is suitable for end-to-end anastomosis, end-to-side anastomosis and lateral anastomosis in digestive tract reconstruction such as esophagus, stomach and intestine.



SCAN THIS CODE
TO VISIT OUR WEBSITE

Address:-A-205 & 208,Premium House,Nr Gandhigram,
Railway Station,Ashram Road,Ahmedabad-380006.

☎ info@fourzone.in | sales@fourzone.in