

SKIN STAPLER

# Rhino Skin Stapler™

*Disposable Skin Stapler*



RHINO SKIN STAPLER

## Rhino Skin Stapler™



RHINO SKIN STAPLER, MAGAZINE AND STAPLE DETAIL

The staple magazine is transparent, so the remaining count is visible at a glance rather than discovered when it runs out, and the angled tip keeps the wound edge in view while firing.

Skin staplers are judged on two things: how fast the incision closes and how the scar looks at six weeks. Rhino addresses both — better apposition of the skin edges for the cosmetic result, and a closure time short enough to shorten the anaesthetic itself on long cases.

— Key Differentiators

- 01

**Transparent magazine**

The remaining staple count is visible at any point during closure.
- 02

**Angled tip**

A clear view of the wound edge, in a light instrument.
- 03

**Edge apposition**

Designed to bring skin edges together cleanly for a better cosmetic result
- 04

**Radiation sterilised**

Supplied sterile, with a matching remover for staple extraction.

— Ordering — three-row haemorrhoid stapler

| CODE        | STAPLE HEIGHT | STAPLES |
|-------------|---------------|---------|
| FZHPILES-34 | 3.8 mm        | 34      |

Same indication as Pile F. Confirm the ordering code with FourZone before raising a purchase order, as the two-row and three-row instruments are listed under a shared reference.

# Secure and stable healthcare products

FourZone Healthcare develops and manufactures devices that improve patient outcomes and advance minimally invasive surgery. Everything in this catalogue is made to one standard, whether it is used in the first minute of a case or the last.

## HEAD OFFICE

---

A-205 & 208, Premium House  
Near Gandhigram Railway Station  
Ashram Road, Ahmedabad 380006  
Gujarat, India

[info@fourzone.in](mailto:info@fourzone.in)  
[sales@fourzone.in](mailto:sales@fourzone.in)  
[fourzone.in](http://fourzone.in)



SCAN TO VISIT