

LAPAROSCOPIC RELOAD

Trizone GC Curved Tip™

Endoscopic Articulating Curved Tip Reload With Funnel And Gecko Claw Technology



TRIZONE GC CURVED TIP RELOAD

Trizone GC Curved Tip™



CURVED-TIP RELOAD, JAW DETAIL AND INSTRUMENT

The curved tip goes around tissue rather than at it, which matters in the low pelvis and behind a mobilised stomach. It also gives a second mechanical benefit: the curve itself resists the tissue sliding forward as the jaw closes

Gecko Claw gripping and the six-row graded staple pattern carry over unchanged.

— Key Differentiators

- 01

Curved tip

Easier to work into and around tissue, with less forward slip during firing.
- 02

Gecko Claw retention

Lateral and longitudinal grip through the full stroke.
- 03

Graded compression

Gradual closure for vascular, medium and thick tissue on the same jaw.
- 04

Consistent staple form

Isometric height retainer keeps formation uniform across the line.

— Ordering — Trizone GC Curved Tip reloads

CODE	CARTRIDGE	STAPLE HEIGHTS (MM)	STAPLE LINE	STAPLES
FZHCL-45W	White	2.4/ 2.6 / 3.0	46 mm	66
FZHCL-45P	Purple	3.4/ 3.6 / 4.0	46 mm	66
FZHCL-45B	Black	4.4/ 4.6 / 4.8	46 mm	66
FZHCL-60W	White	2.4/ 2.6 / 3.0	61 mm	90
FZHCL-60P	Purple	3.4 /3.6 / 4.0	61 mm	90
FZHCL-60B	Black	4.4 / 4.6 / 4.8	61 mm	90

Secure and stable healthcare products

FourZone Healthcare develops and manufactures devices that improve patient outcomes and advance minimally invasive surgery. Everything in this catalogue is made to one standard, whether it is used in the first minute of a case or the last.

HEAD OFFICE

A-205 & 208, Premium House
Near Gandhigram Railway Station
Ashram Road, Ahmedabad 380006
Gujarat, India

info@fourzone.in
sales@fourzone.in
fourzone.in



SCAN TO VISIT